

## PRRT1 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab06832  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information

|              |  |
|--------------|--|
| Immunogen:   | proline-rich transmembrane protein 1   |
| Synonyms:    | Proline-rich transmembrane protein 1 Dispanin subfamily D member 1 (DSPD1) Synapse differentiation-induced protein 4 (SynDIG4) PRRT1 C6orf31 NG5 |
| Observed MW: | 35-40 kDa  |
| Uniprot ID : | Q99946   |

### Application

|                       |                                  |
|-----------------------|----------------------------------|
| Reactivity:           | Human, Mouse, Rat                |
| Tested Application:   | ELISA, WB, IF                    |
| Recommended dilution: | WB: 1:500-1:2000; IF: 1:20-1:200 |
| Image:                |                                  |



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06832(PRRT1 antibody) at dilution of 1:500