

PRDX3 antibody

Product Information

Catalog No.:	FNab06838
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	peroxiredoxin 3
Synonyms:	Thioredoxin-dependent peroxide reductase, mitochondrial Antioxidant protein 1 (AOP-1) HBC189 Peroxiredoxin III (Prx-III) Peroxiredoxin-3 Protein MER5 homolog Thioredoxin-dependent peroxiredoxin 3 PRDX3 AOP1
Observed MW:	28 kDa
Uniprot ID :	P30048

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab06838 (PRX3 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06838 (PRX3 Antibody) at dilution of 1:1000