

## PRDX5 antibody

### Product Information

Catalog No.:	FNab06839
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	peroxiredoxin 5
Synonyms:	Peroxiredoxin-5, mitochondrial Alu corepressor 1 Antioxidant enzyme B166 (AOEB166) Liver tissue 2D-page spot 71B PLP Peroxiredoxin V (Prx-V) Peroxisomal antioxidant enzyme TPx type VI Thioredoxin peroxidase PMP20 Thioredoxin-dependent peroxiredoxin 5 PRDX5 ACR1
Observed MW:	17 kDa
Uniprot ID :	P30044

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab06839(Prx5 antibody) at dilution of 1:50



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab06839(Prx5 antibody) at dilution of 1:600