

PSEN2-Specific antibody

Product Information

Catalog No.:	FNab06853
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	presenilin 2(Alzheimer disease 4)
Synonyms:	Presenilin-2 (PS-2) AD3LP AD5 E5-1 STM-2 Presenilin-2 NTF subunit Presenilin-2 CTF subunit PSEN2 AD4 PS2 PSNL2 STM2
Observed MW:	None
Uniprot ID :	P49810

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab06853(PSEN2-Specific Antibody) at dilution of 1:200