

## TAP2 antibody

### Product Information

Catalog No.:	FNab06854
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	GINS complex subunit 2(Psf2 homolog)
Synonyms:	Antigen peptide transporter 2 (APT2) ATP-binding cassette sub-family B member 3 Peptide supply factor 2 Peptide transporter PSF2 (PSF-2) Peptide transporter TAP2 Peptide transporter involved in antigen processing 2 Really interesting new gene 11 protein (RING11) TAP2 ABCB3 PSF2 RING11 Y1
Observed MW:	25 kDa
Uniprot ID :	Q03519

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab06854(PSF2 antibody) at dilution of 1:2000