

PSMB3 antibody

Product Information

Catalog No.:	FNab06872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The proteasome is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg, Phe, Tyr, Leu, and Glu adjacent to the leaving group at neutral or slightly basic pH. The proteasome has an ATP-dependent proteolytic activity.

Immunogen information


Immunogen:	proteasome(prosome, macropain) subunit, beta type, 3
Synonyms:	Proteasome subunit beta type-3 Proteasome chain 13 Proteasome component C10-II Proteasome theta chain PSMB3
Observed MW:	23-27 kDa
Uniprot ID :	P49720

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab06872(PSMB3 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06872(PSMB3 antibody) at dilution of 1:200