

PSMD4 antibody

Product Information

Catalog No.:	FNab06889
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds and presumably selects ubiquitin-conjugates for destruction. Displays selectivity for longer polyubiquitin chains. Modulates intestinal fluid secretion.

Immunogen information


Immunogen:	proteasome(prosome, macropain) 26S subunit, non-ATPase, 4
Synonyms:	26S proteasome non-ATPase regulatory subunit 4 26S proteasome regulatory subunit RPN10 26S proteasome regulatory subunit S5A Antisecretory factor 1 (AF, ASF) Multiubiquitin chain-binding protein PSMD4 MCB1
Observed MW:	50-55 kDa
Uniprot ID :	P55036

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab06889(PSMD4 Antibody) at dilution of 1:200



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06889(PSMD4 antibody) at dilution of 1:3000