

PSMG3 antibody

Product Information

Catalog No.:	FNab06900
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Chaperone protein which promotes assembly of the 20S proteasome. May cooperate with PSMG1-PSMG2 heterodimers to orchestrate the correct assembly of proteasomes.

Immunogen information

Immunogen:	proteasome(prosome, macropain) assembly chaperone 3
Synonyms:	Proteasome assembly chaperone 3 (PAC-3, hPAC3) PSMG3 C7orf48 PAC3
Observed MW:	None
Uniprot ID :	Q9BT73

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab06900(PSMG3 antibody) at dilution of 1:50