

## PSTPIP2 antibody

### Product Information

Catalog No.:	FNab06904
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds to F-actin. May be involved in regulation of the actin cytoskeleton(By similarity).

### Immunogen information

Immunogen:	proline-serine-threonine phosphatase interacting protein 2
Synonyms:	Proline-serine-threonine phosphatase-interacting protein 2 (PEST phosphatase-interacting protein 2) PSTPIP2
Observed MW:	36 kDa
Uniprot ID :	Q9H939

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with FNab06904(PSTPIP2 antibody) at dilution of 1:500