

PTCD3 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab06909 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |


Background

Immunogen information


| | |
|--------------|--|
| Immunogen: | Pentatricopeptide repeat domain 3 |
| Synonyms: | Small ribosomal subunit protein mS39 28S ribosomal protein S39, mitochondrial (MRP-S39) Pentatricopeptide repeat domain-containing protein 3, mitochondrial Transformation-related gene 15 protein (TRG-15) PTCD3 MRPS39 |
| Observed MW: | 79 kDa |
| Uniprot ID : | Q96EY7 |

Application

| | |
|-----------------------|---------------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:500 - 1:2000; IHC: 1:50 - 1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human stomach using FNab06909(PTCD3 Antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab06909(PTCD3 antibody) at dilution of 1:1000