

PTCH2 antibody

Product Information

Catalog No.:	FNab06911
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May have a role in epidermal development. May act as a receptor for Sonic hedgehog(SHH).

Immunogen information


Immunogen:	patched homolog 2(Drosophila)
Synonyms:	Protein patched homolog 2 (PTC2) PTCH2
Observed MW:	131-140 kDa
Uniprot ID :	Q9Y6C5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06911(PTCH2 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06911(PTCH2 antibody) at dilution of 1:500