

PTHLH antibody

Product Information

Catalog No.:	FNab06922
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	parathyroid hormone-like hormone
Synonyms:	Parathyroid hormone-related protein (PTH-rP, PTHrP) Parathyroid hormone-like protein (PLP) PTHrP[1-36] PTHrP[38-94] Osteostatin Alternative names: PTHrP[107-139] PTHLH PTHRP
Observed MW:	19-24 kDa
Uniprot ID :	P12272

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



various lysates were subjected to SDS PAGE followed by western blot with FNab06922(PTHLH antibody) at dilution of 1:1000