

## PTP4A1 antibody

### Product Information

Catalog No.:	FNab06928
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	protein tyrosine phosphatase type IVA, member 1
Synonyms:	Protein tyrosine phosphatase type IVA 1 PTP(CAAXI) Protein-tyrosine phosphatase 4a1 Protein-tyrosine phosphatase of regenerating liver 1 (PRL-1) PTP4A1 PRL1 PTPCAAX1
Observed MW:	50 kDa
Uniprot ID :	Q93096

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab06928(PTP4A1 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab06928(PTP4A1 antibody) at dilution of 1:800