

PTPN7 antibody

Product Information

Catalog No.:	FNab06940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Protein phosphatase that acts preferentially on tyrosine-phosphorylated MAPK1. Plays a role in the regulation of T and B-lymphocyte development and signal transduction.

Immunogen information


Immunogen:	protein tyrosine phosphatase, non-receptor type 7
Synonyms:	Tyrosine-protein phosphatase non-receptor type 7 Hematopoietic protein-tyrosine phosphatase (HEPTP) Protein-tyrosine phosphatase LC-PTP PTPN7
Observed MW:	38-40 kDa
Uniprot ID :	P35236

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab06940(PTPN7 antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06940(PTPN7 antibody) at dilution of 1:2000