

PTPN9 antibody

Product Information

Catalog No.:	FNab06941
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Protein-tyrosine phosphatase that could participate in the transfer of hydrophobic ligands or in functions of the Golgi apparatus.

Immunogen information

Immunogen:	protein tyrosine phosphatase, non-receptor type 9
Synonyms:	Tyrosine-protein phosphatase non-receptor type 9 Protein-tyrosine phosphatase MEG2 (PTPase MEG2) PTPN9
Observed MW:	68 kDa
Uniprot ID :	P43378

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab06941(PTPN9 Antibody) at dilution of 1:1000