

PTPRO antibody

Product Information

Catalog No.:	FNab06946
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Possesses tyrosine phosphatase activity. Plays a role in regulating the glomerular pressure/filtration rate relationship through an effect on podocyte structure and function(By similarity).

Immunogen information

Immunogen:	protein tyrosine phosphatase, receptor type, O
Synonyms:	Receptor-type tyrosine-protein phosphatase O (R-PTP-O) Glomerular epithelial protein 1 Protein tyrosine phosphatase U2 (PTP-U2, PTPase U2) PTPRO GLEPP1 PTPU2
Observed MW:	None
Uniprot ID :	Q16827

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA