

PUS3 antibody

Product Information

Catalog No.:	FNab06959
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Formation of pseudouridine at position 39 in the anticodon stem and loop of transfer RNAs.

Immunogen information

Immunogen:	pseudouridylate synthase 3
Synonyms:	tRNA pseudouridine(38/39) synthase tRNA pseudouridine synthase 3 tRNA pseudouridylate synthase 3 tRNA-uridine isomerase 3 PUS3
Observed MW:	50-52 kDa
Uniprot ID :	Q9BZE2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab06959(PUS3 antibody) at dilution of 1:500