

NECTIN4 antibody

Product Information

Catalog No.:	FNab06962
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	poliovirus receptor-related 4
Synonyms:	Nectin-4 Ig superfamily receptor LNIR Nectin cell adhesion molecule 4 Poliovirus receptor-related protein 4 Processed poliovirus receptor-related protein 4 NECTIN4 LNIR PRR4 PVRL4
Observed MW:	43 kDa, 66 kDa
Uniprot ID :	Q96NY8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200; IHC: 1:200-1:2000
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab06962(PVRL4 antibody) at dilution of 1:1000



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06962(PVRL4 antibody) at dilution of 1:2000