

PWP2 antibody

Product Information

Catalog No.:	FNab06964
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	PWP2 periodic tryptophan protein homolog(yeast)
Synonyms:	Periodic tryptophan protein 2 homolog PWP2 PWP2H
Observed MW:	102 kDa
Uniprot ID :	Q15269

Application

Reactivity:	Human
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06964(PWP2 Antibody) at dilution of 1:600