

## NR1I2 antibody

### Product Information


Catalog No.:	FNab06967
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	nuclear receptor subfamily 1, group I, member 2
Synonyms:	Nuclear receptor subfamily 1 group I member 2 Orphan nuclear receptor PAR1 Orphan nuclear receptor PXR Pregnane X receptor Steroid and xenobiotic receptor (SXR) NR1I2 PXR
Observed MW:	53 kDa
Uniprot ID :	O75469

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse testis using FNab06967(PXR antibody) at dilution of 1:100



Rat liver tissue were subjected to SDS PAGE followed by western blot with FNab06967(PXR antibody) at dilution of 1:1000