

PTK2B antibody

Product Information

Catalog No.:	FNab06980
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	PTK2B protein tyrosine kinase 2 beta
Synonyms:	Protein-tyrosine kinase 2-beta Calcium-dependent tyrosine kinase (CADTK) Calcium-regulated non-receptor proline-rich tyrosine kinase Cell adhesion kinase beta (CAK-beta, CAKB) Focal adhesion kinase 2 (FADK 2) Proline-rich tyrosine kinase 2 Related adhesion focal tyrosine kinase (RAFTK) PTK2B FAK2 PYK2 RAFTK
Observed MW:	112-115 kDa
Uniprot ID :	Q14289

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human medulloblastoma using FNab06980(PTK2B antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06980(PTK2B antibody) at dilution of 1:1000