

QDPR antibody

Product Information

Catalog No.:	FNab06982
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The product of this enzyme, tetrahydrobiopterin(BH-4), is an essential cofactor for phenylalanine, tyrosine, and tryptophan hydroxylases.

Immunogen information


Immunogen:	quinoid dihydropteridine reductase
Synonyms:	Dihydropteridine reductase HDHPR Quinoid dihydropteridine reductase Short chain dehydrogenase/reductase family 33C member 1 QDPR DHPR SDR33C1
Observed MW:	28-30 kDa
Uniprot ID :	P09417

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab06982(QDPR antibody) at dilution of 1:100



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab06982(QDPR antibody) at dilution of 1:1000