

## QPRT antibody

### Product Information

Catalog No.:	FNab06984
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the catabolism of quinolinic acid(QA).

### Immunogen information

Immunogen:	quinolinate phosphoribosyltransferase
Synonyms:	Nicotinate-nucleotide pyrophosphorylase [carboxylating] Quinolinate phosphoribosyltransferase [decarboxylating] (QAPRTase, QPRTase) QPRT
Observed MW:	31 kDa
Uniprot ID :	Q15274

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06984( QPRT Antibody) at dilution of 1:1000