

RAB11FIP5 antibody

Product Information

Catalog No.:	FNab06993
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Rab effector involved in protein trafficking from apical recycling endosomes to the apical plasma membrane. Involved in insulin granule exocytosis. May regulate V-ATPase intracellular transport in response to extracellular acidosis.

Immunogen information

Immunogen:	RAB11 family interacting protein 5(class I)
Synonyms:	Rab11 family-interacting protein 5 (Rab11-FIP5) Gamma-SNAP-associated factor 1 (Gaf-1) Phosphoprotein pp75 Rab11-interacting protein Rip11 RAB11FIP5 GAF1 KIAA0857 RIP11
Observed MW:	66-72 kDa
Uniprot ID :	Q9BXF6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab06993(RAB11FIP5 antibody) at dilution of 1:500