

RAB11-Specific antibody

Product Information


Catalog No.:	FNab06994
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	RAB11A, member RAS oncogene family
Synonyms:	Ras-related protein Rab-11A (Rab-11) YL8 RAB11A RAB11
Observed MW:	24 kDa
Uniprot ID :	P62491

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-RAB11A-Specific
(IP:FNab06994, 3ug; Detection:FNab06994 1:700)
with SH-SY5Y cells lysate 3000ug.



MCF-7 cells were subjected to SDS PAGE
followed by western blot with 20229-1-AP
(RAB11A-Specific antibody) at dilution of 1:1000