

RAB33A antibody

Product Information


Catalog No.:	FNab07017
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	RAB33A, member RAS oncogene family
Synonyms:	Ras-related protein Rab-33A Small GTP-binding protein S10 RAB33A RABS10
Observed MW:	27 kDa
Uniprot ID :	Q14088

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab07017(RAB33A antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07017(RAB33A antibody) at dilution of 1:600