

## RAB34 antibody

### Product Information

Catalog No.:	FNab07018
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Protein transport. Involved in the redistribution of lysosomes to the peri-Golgi region(By similarity). Plays a role in the maturation of phagosomes that engulf pathogens, such as S.aureus and M.tuberculosis. Plays a role in the fusion of phagosomes with lysosomes.

### Immunogen information


Immunogen:	RAB34, member RAS oncogene family
Synonyms:	Ras-related protein Rab-34 Ras-related protein Rab-39 Ras-related protein Rah RAB34 RAB39 RAH
Observed MW:	29 kDa
Uniprot ID :	Q9BZG1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IP: 1:200-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab07018(RAB39 Antibody) at dilution of 1:200. heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07018(RAB39 Antibody) at dilution of 1:2000