

RAB37 antibody

Product Information


Catalog No.:	FNab07020
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	RAB37, member RAS oncogene family
Synonyms:	Ras-related protein Rab-37 RAB37
Observed MW:	25 kDa
Uniprot ID :	Q96AX2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab07020(RAB37 antibody) at dilution of 1:50



Sp2/0 cells were subjected to SDS PAGE followed by western blot with FNab07020(RAB37 antibody) at dilution of 1:100