

RAB3B antibody

Product Information

Catalog No.:	FNab07026
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Protein transport. Probably involved in vesicular traffic(By similarity).

Immunogen information


Immunogen:	RAB3B, member RAS oncogene family
Synonyms:	Ras-related protein Rab-3B RAB3B
Observed MW:	25 kDa
Uniprot ID :	P20337

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab07026(RAB3B antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07026(RAB3B antibody) at dilution of 1:800