

## RAB3C-specific antibody

### Product Information


Catalog No.:	FNab07028
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	RAB3C, member RAS oncogene family
Synonyms:	Ras-related protein Rab-3C RAB3C
Observed MW:	26 kDa
Uniprot ID :	Q96E17

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:1000; IF: 1:10-1:100
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07028(RAB3C-specific antibody) at dilution of 1:300