

RAB5A-Specific antibody

Product Information

Catalog No.:	FNab07039
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. The antibody is specific to RAB5A. It has no cross reaction to RAB5B or RAB5C.

Immunogen information


Immunogen:	RAB5A, member RAS oncogene family
Synonyms:	Ras-related protein Rab-5A RAB5A RAB5
Observed MW:	24 kDa
Uniprot ID :	P20339

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab07039(RAB5A-Specific Antibody) at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07039(RAB5A-Specific antibody) at dilution of 1:500