

RAB5B antibody

Product Information

Catalog No.:	FNab07040
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Protein transport. Probably involved in vesicular traffic(By similarity).

Immunogen information

Immunogen:	RAB5B, member RAS oncogene family
Synonyms:	Ras-related protein Rab-5B RAB5B
Observed MW:	27 kDa
Uniprot ID :	P61020

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab07040(RAB5B antibody) at dilution of 1:1000