

## RAB9A antibody

### Product Information

Catalog No.:	FNab07048
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the transport of proteins between the endosomes and the trans Golgi network. Involved in the recruitment of SGSM2 to melanosomes and is required for the proper trafficking of melanogenic enzymes TYR, TYRP1 and DCT/TYRP2 to melanosomes in melanocytes.

### Immunogen information


Immunogen:	RAB9A, member RAS oncogene family
Synonyms:	Ras-related protein Rab-9A RAB9A RAB9
Observed MW:	23 kDa
Uniprot ID :	P51151

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab07048(RAB9A antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07048(RAB9A antibody) at dilution of 1:300