

RAB9B antibody

Product Information

Catalog No.:	FNab07050
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the transport of proteins between the endosomes and the trans Golgi network.

Immunogen information

Immunogen:	RAB9B, member RAS oncogene family
Synonyms:	Ras-related protein Rab-9B Rab-9-like protein (Rab-9L) RAB9B RAB9L
Observed MW:	23-25 kDa
Uniprot ID :	Q9NP90

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07050(RAB9B antibody) at dilution of 1:300