

## RABL2B antibody

### Product Information


Catalog No.:	FNab07062
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	RAB, member of RAS oncogene family-like 2B
Synonyms:	Rab-like protein 2B RABL2B
Observed MW:	26 kDa
Uniprot ID :	Q9UNT1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab07062(RABL2B antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07062(RABL2B antibody) at dilution of 1:500