

RALBP1 antibody

Product Information

Catalog No.:	FNab07095
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

RALBP1 plays a role in receptor-mediated endocytosis and is a downstream effector of the small GTP-binding protein RAL (see RALA; MIM 179550). Small G proteins, such as RAL, have GDP-bound inactive and GTP-bound active forms, which shift from the inactive to the active state through the action of RALGDS (MIM 601619), which in turn is activated by RAS (see HRAS; MIM 190020) (summary by Feig, 2003 [PubMed 12888294]).


Immunogen information

Immunogen:	ralA binding protein 1
Synonyms:	RalA-binding protein 1 (RalBP1) 76 kDa Ral-interacting protein Dinitrophenyl S-glutathione ATPase (DNP-SG ATPase) Ral-interacting protein 1 RALBP1 RLIP RLIP1 RLIP76
Observed MW:	100 kDa
Uniprot ID :	Q15311


Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Image:



Immunohistochemistry of paraffin-embedded rat heart using FNab07095(RALBP1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07095(RALBP1 antibody) at dilution of 1:1000