

## RAP1GAP antibody

### Product Information

Catalog No.:	FNab07110
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


GTPase activator for the nuclear Ras-related regulatory protein RAP-1A(KREV-1), converting it to the putatively inactive GDP-bound state.

### Immunogen information


Immunogen:	RAP1 GTPase activating protein
Synonyms:	Rap1 GTPase-activating protein 1 (Rap1GAP, Rap1GAP1) RAP1GAP KIAA0474 RAP1GA1
Observed MW:	85-95 kDa
Uniprot ID :	P47736

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab07110(RAP1GAP Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07110(RAP1GAP antibody) at dilution of 1:2000