

RAPGEF5 antibody

Product Information

Catalog No.:	FNab07115
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Guanine nucleotide exchange factor(GEF) for RAP1A, RAP2A and MRAS/M-Ras-GTP. Its association with MRAS inhibits Rap1 activation.

Immunogen information

Immunogen:	Rap guanine nucleotide exchange factor(GEF) 5
Synonyms:	Rap guanine nucleotide exchange factor 5 Guanine nucleotide exchange factor for Rap1 M-Ras-regulated Rap GEF (MR-GEF) Related to Epac (Repac) RAPGEF5 GFR KIAA0277 MRGEF
Observed MW:	68-75 kDa
Uniprot ID :	Q92565

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07115(RAPGEF5 antibody) at dilution of 1:500