

## RAPGEFL1 antibody

### Product Information

Catalog No.:	FNab07116
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probable guanine nucleotide exchange factor(GEF).

### Immunogen information

Immunogen:	Rap guanine nucleotide exchange factor(GEF)-like 1
Synonyms:	Rap guanine nucleotide exchange factor-like 1 Link guanine nucleotide exchange factor II (Link GEFII) RAPGEFL1
Observed MW:	73 kDa, 59 kDa
Uniprot ID :	Q9UHV5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07116(RAPGEFL1 antibody) at dilution of 1:500