

RASSF4 antibody

Product Information

Catalog No.:	FNab07136
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Potential tumor suppressor. May act as a KRAS effector protein. May promote apoptosis and cell cycle arrest.

Immunogen information


Immunogen:	Ras association(RalGDS/AF-6) domain family member 4
Synonyms:	Ras association domain-containing protein 4 RASSF4
Observed MW:	37 kDa
Uniprot ID :	Q9H2L5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab07136(RASSF4 antibody) at dilution of 1:100



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07136(RASSF4 Antibody) at dilution of 1:1000