

RASSF7 antibody

Product Information

Catalog No.:	FNab07138
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Ras association(RalGDS/AF-6) domain family(N-terminal) member 7
Synonyms:	Ras association domain-containing protein 7 HRAS1-related cluster protein 1 RASSF7 C11orf13 HRC1
Observed MW:	40 kDa
Uniprot ID :	Q02833

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-12 cells were subjected to SDS PAGE followed by western blot with FNab07138(RASSF7 antibody) at dilution of 1:500