

RBBP9 antibody

Product Information

Catalog No.:	FNab07156
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7E10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in the transformation process due to its capacity to confer resistance to the growth-inhibitory effects of TGF-beta1 through interaction with retinoblastoma and the subsequent displacement of E2F-1.

Immunogen information


Immunogen:	retinoblastoma binding protein 9
Synonyms:	Serine hydrolase RBBP9 B5T-overexpressed gene protein (Protein BOG) Retinoblastoma-binding protein 10 (RBBP-10) Retinoblastoma-binding protein 9 (RBBP-9) RBBP9 BOG RBBP10
Observed MW:	22 kDa
Uniprot ID :	O75884

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab07156(RBBP9 antibody) at dilution of 1:50



A549 cells were subjected to SDS PAGE followed by western blot with FNab07156(RBBP9 antibody) at dilution of 1:500