

RBM10 antibody

Product Information

Catalog No.:	FNab07158
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in post-transcriptional processing, most probably in mRNA splicing. Binds to RNA homopolymers, with a preference for poly(G) and poly(U) and little for poly(A)(By similarity).

Immunogen information

Immunogen:	RNA binding motif protein 10
Synonyms:	RNA-binding protein 10 G patch domain-containing protein 9 RNA-binding motif protein 10 RNA-binding protein S1-1 (S1-1) RBM10 DXS8237E GPATC9 GPATCH9 KIAA0122
Observed MW:	115 kDa
Uniprot ID :	P98175

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07158(RBM10 antibody) at dilution of 1:1000