

RBM15B antibody

Product Information

Catalog No.:	FNab07164
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function in the regulation of alternative or illicit splicing.

Immunogen information

Immunogen:	RNA binding motif protein 15B
Synonyms:	Putative RNA-binding protein 15B One-twenty two protein 3 (HsOTT3, HuOTT3) RNA-binding motif protein 15B RBM15B OTT3
Observed MW:	100-120 kDa, 70 kDa
Uniprot ID :	Q8NDT2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07164(RBM15B antibody) at dilution of 1:2000