

SCAF8 antibody

Product Information

Catalog No.:	FNab07165
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	RNA binding motif protein 16
Synonyms:	SR-related and CTD-associated factor 8 CDC5L complex-associated protein 7 RNA-binding motif protein 16 SCAF8 CCAP7 KIAA1116 RBM16
Observed MW:	170 kDa
Uniprot ID :	Q9UPN6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:1000-1:10000; IP:1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab07165(RBM16 antibody) at dilution of 1:100



HT-1080 cells were subjected to SDS PAGE followed by western blot with FNab07165(RBM16 antibody) at dilution of 1:2000