

RBM42 antibody

Product Information

Catalog No.:	FNab07179
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds(via the RRM domain) to the 3'-untranslated region(UTR) of CDKN1A mRNA.

Immunogen information

Immunogen:	RNA binding motif protein 42
Synonyms:	RNA-binding protein 42 RNA-binding motif protein 42 RBM42
Observed MW:	65-70 kDa
Uniprot ID :	Q9BTD8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07179(RBM42 antibody) at dilution of 1:500