

RBP1 antibody

Product Information

Catalog No.:	FNab07188
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Intracellular transport of retinol.

Immunogen information


Immunogen:	retinol binding protein 1, cellular
Synonyms:	Retinol-binding protein 1 Cellular retinol-binding protein (CRBP) Cellular retinol-binding protein I (CRBP-I) RBP1 CRBP1
Observed MW:	14-21 kDa
Uniprot ID :	P09455

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNa07188(RBP1 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNa07188(RBP1 antibody) at dilution of 1:1000