

RBPJL antibody

Product Information

Catalog No.:	FNab07197
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Putative transcription factor, which cooperates with EBNA2 to activate transcription.

Immunogen information

Immunogen:	recombination signal binding protein for immunoglobulin kappa J region-like
Synonyms:	Recombining binding protein suppressor of hairless-like protein Transcription factor RBP-L RBPJL RBPL RBPSUHL
Observed MW:	57-69 kDa
Uniprot ID :	Q9UBG7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07197(RBPJL Antibody) at dilution of 1:1500