

RCAN1 antibody

Product Information

Catalog No.:	FNab07201
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.

Immunogen information


Immunogen:	regulator of calcineurin 1
Synonyms:	Calcipressin-1 Adapt78 Down syndrome critical region protein 1 Myocyte-enriched calcineurin-interacting protein 1 (MCIP1) Regulator of calcineurin 1 RCAN1 ADAPT78 CSP1 DSC1 DSCR1
Observed MW:	25-29 kDa
Uniprot ID :	P53805

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab07201(RCAN1 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07201(RCAN1 antibody) at dilution of 1:500