

EBAG9 antibody

Product Information

Catalog No.:	FNab07204
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9E2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	estrogen receptor binding site associated, antigen, 9
Synonyms:	Receptor-binding cancer antigen expressed on SiSo cells Cancer-associated surface antigen RCAS1 Estrogen receptor-binding fragment-associated gene 9 protein EBAG9 RCAS1
Observed MW:	28 kDa
Uniprot ID :	O00559

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab07204(EBAG9 Antibody) at dilution of 1:50 Heat mediated antigen retrieved with Citric acid buffer, pH6.0.



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07204(EBAG9 Antibody) at dilution of 1:1000